

Branched Chain Amino Acid (Leu, Ile, Val) Colorimetric Assay Kit

(Catalog #K564-100; 100 Reactions; Store kit at -20°C)

I. Introduction:

The branched-chain amino acids or BCAA's, refer to the amino acids with non-linear aliphatic side-chains, namely leucine, isoleucine and valine. These three essential amino acids make up approximately 1/3 of skeletal muscle in the human body. BCAA's are currently used clinically to aid in the recovery of burn victims, as well as for strength supplementation for athletes. BCAA's, primarily Leu, can stimulate insulin secretion. The BCAA's have also been implicated in a wide range of other physiological effects. BioVision's BCAA Assay Kit provides a simple convenient means of measuring the BCAA's in a variety of biological samples. The kit utilizes an enzyme assay in which BCAA is oxidatively deaminated, producing NADH which reduces the probe, generating a colored product (λ_{max} = 450 nm). BioVision's BCAA kit measures BCAA's in the range of 0 to 10 nmol per sample with a detection limit of ~0.2 nmol (~10 μ M BCAA in sample). BCAA's are present in serum ~ 0.1-0.4 mM each (~0.125-1.5 mM combined).

II. Kit Contents:

Components	K564-100	Cap Code	Part No.
BCAA Assay Buffer	25 ml	WM	K564-100-1
BCAA Enzyme Mix	lyophilized	Green	K564-100-2
WST Substrate Mix	lyophilized	Red	K564-100-3
Leu Standard (1 μ mol)	100 μ l	Yellow	K564-100-4

III. Storage and Handling:

Store the kit at -20°C, protect from light. Allow Assay Buffer to warm to room temperature before use. Briefly centrifuge vials prior to opening. Read the entire protocol before performing the assay.

IV. Reagent Reconstitution and General Consideration:

BCAA Enzyme Mix: Dissolve with 220 μ l BCAA Assay Buffer. Pipette up and down to dissolve. Stable at 4°C for two months.

WST Substrate Mix: Dissolve with 220 μ l of dH₂O before use. Mix well, store at 4°C, protect from light. Stable for 2 months.

Leucine Standard: Ready to use as supplied. Store at 4°C.

V. BCAA Assay Protocol:

- Standard Curve:** Dilute 10 μ l of the 10 mM Leucine Standard with 90 μ l dH₂O to generate 1 mM Leucine standard. Add 0, 2, 4, 6, 8, 10 μ l of the diluted Standard into a 96-well plate to generate 0, 2, 4, 6, 8, 10 nmol/well standard. Bring the volume to 50 μ l with Assay Buffer.

- Sample Preparation:**

Tissue (20 mg) or cells (2×10^6) can be homogenized with 100 μ l Assay buffer. Centrifuge at 15,000g for 10 minutes to remove cell debris and other insoluble materials. Add samples to sample wells in a 96-well plate and bring the volume to 50 μ l/well with Assay Buffer. We suggest testing several doses of your sample to make sure the readings are within the standard curve range. Typical volume for serum samples should be in the range of 1 – 20 μ l.

3. Reaction Mix:

Mix enough reagents for the number of assays to be performed. For each well, prepare a total 50 μ l Reaction Mix containing:

Amino Acid Measurement		Bkgd Control
Assay Buffer	46 μ l	48 μ l
Enzyme Mix	2 μ l	-----
WST Substrate Mix	2 μ l	2 μ l

Add 50 μ l of the Reaction Mix to each well containing the leucine standard and test samples. Mix well. Incubate the reaction for 30 min at room temperature, protect from light. NADH and NADPH can generate significant background. If these compounds are suspected of being in your sample at significant concentration, perform a simple background control by replacing the Enzyme Mix with 2 μ l Assay Buffer. The background reading should be subtracted from the BCAA test sample readings.

- Measure O.D. at 450 nm in a microplate reader


- Calculation:** Correct background by subtracting the value derived from the 0 BCAA standards from all readings (The background reading can be significant and must be subtracted from sample readings). Plot standard curve. Apply sample readings to the standard curve. BCAA concentrations of the test samples can then be calculated:

$$C = S_a/S_v \text{ (nmol/}\mu\text{l, or mM)}$$

Where:

S_a = BCAA content of unknown samples (nmol) from standard curve,
 S_v = sample volume (μ l) added into the assay wells.

BCAA molecular weights are: Leu 131.18, Ile 131.18, Val 117.15 g/mol.



Leucine Assay performed according to this protocol

RELATED PRODUCTS:

NAD(P)/NAD(P)H Quantification Kit	ADP/ATP Ratio Assay Kit
Ascorbic Acid Quantification Kit	Glutathione Detection Kit
Total Antioxidant Capacity (TAC) Assay Kit	Fatty Acid Assay Kit
Ethanol Assay Kit	Uric Acid Assay Kit
Pyruvate Assay Kit	Lactate Assay Kit I & II
Creatinine Assay Kit	Nitric Oxide Assay Kit
Ammonia Assay Kit	Free Glycerol Assay Kit
Triglyceride Assay Kit	Hemin Assay Kit
Alanine Assay Kit	Glucose Assay Kit
Sarcosine Assay Kit	L-Amino Acid Assay Kit
Phenylalanine Assay Kit	Cholesterol Assay Kit

FOR RESEARCH USE ONLY! Not to be used on humans.

GENERAL TROUBLESHOOTING GUIDE:

Problems	Cause	Solution
Assay not working	• Use of ice-cold assay buffer • Omission of a step in the protocol • Plate read at incorrect wavelength • Use of a different 96-well plate	• Assay buffer must be at room temperature • Refer and follow the data sheet precisely • Check the wavelength in the data sheet and the filter settings of the instrument • Fluorescence: Black plates ; Luminescence: White plates ; Colorimeters: Clear plates
Samples with erratic readings	• Use of an incompatible sample type • Samples prepared in a different buffer • Samples were not deproteinized (if indicated in datasheet) • Cell/ tissue samples were not completely homogenized • Samples used after multiple free-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples	• Refer data sheet for details about incompatible samples • Use the assay buffer provided in the kit or refer data sheet for instructions • Use the 10 kDa spin cut-off filter or PCA precipitation as indicated • Use Dounce homogenizer (increase the number of strokes); observe for lysis under microscope • Aliquot and freeze samples if needed to use multiple times • Troubleshoot if needed, deproteinize samples • Use fresh samples or store at correct temperatures till use
Lower/ Higher readings in Samples and Standards	• Improperly thawed components • Use of expired kit or improperly stored reagents • Allowing the reagents to sit for extended times on ice • Incorrect incubation times or temperatures • Incorrect volumes used	• Thaw all components completely and mix gently before use • Always check the expiry date and store the components appropriately • Always thaw and prepare fresh reaction mix before use • Refer datasheet & verify correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly
Readings do not follow a linear pattern for Standard curve	• Use of partially thawed components • Pipetting errors in the standard • Pipetting errors in the reaction mix • Air bubbles formed in well • Standard stock is at an incorrect concentration • Calculation errors • Substituting reagents from older kits/ lots	• Thaw and resuspend all components before preparing the reaction mix • Avoid pipetting small volumes • Prepare a master reaction mix whenever possible • Pipette gently against the wall of the tubes • Always refer the dilutions in the data sheet • Recheck calculations after referring the data sheet • Use fresh components from the same kit
Unanticipated results	• Measured at incorrect wavelength • Samples contain interfering substances • Use of incompatible sample type • Sample readings above/below the linear range	• Check the equipment and the filter setting • Troubleshoot if it interferes with the kit • Refer data sheet to check if sample is compatible with the kit or optimization is needed • Concentrate/ Dilute sample so as to be in the linear range
Note: The most probable list of causes is under each problem section. Causes/ Solutions may overlap with other problems.		